SARS-CoV-2
Antigen

ID: _____

C
T

S

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INTENDED USE

This rapid test kit is intended for the qualitative detection of SARS-CoV-2 infection from patients. It is for professional use only. It is an aid in the diagnosis of the patients with suspected SARS-CoV-2 infection in conjunction with clinical presentation and results of other laboratory tests. Results from this test kit should not be used as the sole basis for diagnosis.

The test provides preliminary test results. Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for treatment or other management decision.

PRINCIPLE

This kit is an immunochromatography assay. According to the gold immunochromatographic test principle, double antibody sandwich method was used to detect SARS-CoV-2 antigen in the samples. When there is virus antigen presence in the sample, the antigen binds with the corresponding colloidal gold monoclonal antibody and the coated monoclonal antibody at the detection line to form a compound and then condenses into a red band, indicating a positive result. If there is no antigen in the sample, complex cannot be formed at the detection line, and no red band is shown, indicating negative result.

Whether the sample contains antigen or not, the gold monoclonal antibody will bind to the enveloped antibody at the quality control line, form a compound and condense into a red band.

KIT COMPONENTS

- 20 Test cassettes
- 2 Sample extraction
- 20 Sample tubes
- 20 single packaged Swabs
- 1 Tube stand
- 1 Instruction for use

STORAGE AND EXPIRY

- Store as packaged in the sealed pouch at 2-30°C, avoid hot and sunshine, dry place, valid for 12 months.
- Do not freeze. Some protective measures should be taken in hot summer and cold winter to avoid high temperature or freeze-thaw.
- Do not open the inner packaging until ready, it must be used in one hour if opened (Humidity≤60%, Temp: 20°C-30°C).
- Please use immediately when the humidity>60%.

SAMPLE PRESERVATION

Sample of human nasopharyngeal swabs and oropharyngeal swabs should be processed as soon as possible after sample collection. If the test cannot be performed immediately, the sample should be stored in a sealed state, stored at 2~8°C for 8 hours, and stored below -20°C for

1 month. Long-term storage is not recommended.

SAMPLE COLLECTION

By nasopharyngeal swab: The sampler gently holds the head of the person to be collected with one hand, holds the swab with the other hand, sticks the swab to the nostril to enter, and slowly penetrates backwards along the bottom of the lower nasal passage, so as not to exert too much force to avoid traumatic hemorrhage. When the tip of the swab reaches the posterior wall of the nasopharynx cavity, gently rotate it (make sure the swab contains both cells and slime), and then slowly take out the swab.

By oropharyngeal swab: The head of the person to be collected is slightly tilted and his mouth is wide open, exposing the pharyngeal tonsils on both sides. Wipe the swab across the root of the tongue.

Wipe the pharyngeal tonsils on both sides of the person to be collected back and forth with a little force for at least 3 times, and then wipe up and down the posterior pharyngeal wall for at least 3 times.

TEST STEPS

Instructions must be read entirely before taking the test. Bring test cassettes and test components at room temperature before performing the test.

1. Add 23 drops sample extraction (550 µl) into the sample tube.

2. Dip the swab after collecting the sample into the sample extraction liquid, fully immerse the tip of the swab, rotate and squeeze the swab 10 times, Roll the swab head against the inner wall of the tube as you remove it. Try to release as much liquid as possible. Dispose the used swab in accordance with biohazard waste disposal protocol.

3. Take out the test cassette from the sealed pouch, place it on a clean and level surface with the sample port well up. Name the test cassette.

4. Apply 2 full drops of the treated sample (60µl-70µl) vertically into the sample well of the test cassette.

5. Observe the test results immediately within 15~20 minutes, the result is invalid over 20 minutes.

Note: When using viral transport medium (VTM), it is important to ensure that the VTM containing the sample is warmed to room temperature. Cold samples will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold sample to room temperature.

INTERPRETATION OF RESULTS

POSITIVE: Two (2) distinct colored lines appear. One line should be in the control region (C) and the other line should be in the test region (T).

NEGATIVE: One (1) colored line appears in the control region (C). No apparent colored line appears in the test region (T). The negative result does not indicate the absence of analytes in the sample.

INVALID: No colored lines appear, or control line fails to appear, indicating that the operator error or reagent failure. Verify the test procedure and repeat the test.

LIMITATIONS

1. This reagent is a qualitative detection reagent, which cannot determine the exact content of antigen.
2. The test results of this reagent are only for the reference of clinicians and should not be taken as the sole basis for clinical diagnosis and treatment. Clinical management of patients should be considered in the light of their symptoms/signs, medical history, other laboratory tests and treatment responses.
3. Restricted by antigen detection reagent method, the lowest detection limit (sensitivity analysis) is generally lower than that of nucleic acid detection, so the researchers deal with negative result to give more attention, should be combined with other test results comprehensive judgment, advice to doubt the negative result of nucleic acid detection or virus isolation culture identification method for review.
4. False negative results may be caused by unreasonable sample collection, transport and treatment, and low viral load in samples.

PERFORMANCE CHARACTERISTICS

1. Sensitivity and Specificity

The verification of naso-oropharyngeal samples is carried out with enterprise reference and the results meet the following requirements:

The National Reference Panel for SARS-Cov-2 Antigen Detection Kit of the National Institutes for Food and Drug Control was used for verification, and the results met the following requirements:

Coincidence rate of negative reference products: National negative reference materials (N1-N20) were tested, and the results were all negative, with a coincidence rate of 100%. Coincidence rate of positive reference products: tested with national positive reference products (P1-P8), the results are all positive, and the coincidence rate is 100%.

Repeatability: The 10 test results of R1 and R2 should be positive, and the color rendering is uniform.

Clinical Sensitivity and Specificity:

		RT-PCR		Total
		Positive	Negative	
SARS-CoV-2 Antigen Rapid Test	Positive	97	6	103
	Negative	3	294	297
Total		100	300	400

Diagnostic sensitivity: 97,0% (91,5% - 99,4%) *

Diagnostic specificity: 98,0% (95,7% - 99,3%) *

Total coincidence rate: 97,8% (95,8% - 99,0%) *

*95 % CI (Confidence Interval)

2. Limit of detection:

LOD concentration	2.4 * 10 ² TCID ₅₀ /mL
-------------------	--

3. Cross-reactivity

With human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, adenovirus, human metapneumovirus, parainfluenza virus 1, parainfluenza virus 2, parainfluenza virus 3, parainfluenza virus 4, influenza A, type B Influenza, enterovirus, respiratory syncytial virus, rhinovirus, SARS coronavirus, MERS coronavirus, Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bordetella pertussis, Mycoplasma pneumoniae, pneumonia Chlamydia, Legionella pneumophila, etc. have no cross reaction (see table in the German IFU).

4. Interfering

Common interfering substances in the sample, such as blood, mucin, and pus, have no effect on the test results.

PRECAUTIONS

1. For IN VITRO diagnostic use only.
2. Reagents should be used as soon as possible after opened. This reagent cannot be reused for disposable.
3. The test device should remain in the sealed pouches until use.
4. All specimens and reagents should be considered potentially hazardous and handled in the same manner as an infectious agent after use.
5. Don't use after the expiration date.

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Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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product group	product for self-testing	art.no.	language (if applicable country group LG)
basic product	NANOREPRO SARS-CoV-2 Antigen Schnelltest	B66000	DE/EN

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***immunochemographic rapid test for professional use for the qualitative detection of SARS-CoV-2
Infection from patients***

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AS, 

Anne Schelling

Head of Quality Management
Safety Officer for Medical Devices
Head of Regulatory Affairs

Performance Evaluation of NANOREPRO SARS-CoV-2 Antigen Rapid Test

1. Clinical Evaluation

1.1 Objective

To determine the clinical sensitivity and specificity of NANOREPRO SARS-CoV-2 Antigen Rapid Test by comparing with RT-PCR for the intended use of detecting SARS-CoV-2 viral antigens in patients within two weeks of symptom onset.

1.2 Parameter

Clinical evaluation was performed to compare the results obtained by NANOREPRO SARS-CoV-2 Antigen Rapid Test and a real-time Polymerase Chain Reaction (RT-PCR) comparator assay. Patients who presented with symptoms within 7 days were included in the study.

400 clinical samples were selected for the study. There were 100 clinical samples from the Individuals who were suspected of COVID-19 within 7 days of symptom and 300 negative samples. The NANOREPRO SARS-CoV-2 Antigen Rapid Test and the PCR test were detected simultaneously, and the positive coincidence rate as well as negative coincidence rate was calculated.

1.3 Result

		RT-PCR		Total
		Positive	Negative	
NANOREPRO SARS-CoV-2 Antigen Rapid Test	Positive	97	6	103
	Negative	3	294	297
	Total	100	300	400

Relative Sensitivity: 97.00%

Relative Specificity:98.00%

1.4 Conclusion

The above results demonstrate that the NANOREPRO SARS-CoV-2 Antigen Rapid Test has relative sensitivity of 97% and relative specificity of 98%

2. Cross Reactivity

2.1 Objective

To determine the cross reactivity of NANOREPRO SARS-CoV-2 Antigen Rapid Test

2.2 Parameter

In order to determine cross reactivity samples which are positive for following organisms were collected and tested with NANOREPRO SARS-CoV-2 Antigen Rapid Test

Human coronavirus (229E/OC43/NL63/SARS/MERS)	Candida albicans
Influenza A	Streptococcus pneumoniae
Influenza B	Streptococcus pyogenes
Parainfluenza (1/2/3/4)	Adenovirus
Rhinovirus	Mycoplasma pneumoniae
Chlamydia pneumoniae	Legionella pneumophila
RSV	Bordetella pertussis
Enterovirus	Pooled human nasal wash
Haemophilus influenza	
Human Metapneumovirus	

2.3 Result / Conclusion

No cross reactivity was observed and the samples mentioned above showed negative result when tested with SARS-CoV-2 Antigen Rapid Test

3. Interference Study

3.1 Objective

To determine interference study of NANOREPRO SARS-CoV-2 Antigen Rapid Test

3.2 Result / Conclusion

The performance of NANOREPRO SARS-CoV-2 Antigen Rapid Test is not interfered by blood, mucin and pus.

4. Analytical Sensitivity (Limit of Detection)

4.1 Objective

This study was performed to determine the analytical sensitivity of NANOREPRO SARS-CoV-2 Antigen Rapid Test.

4.2 Result

The limit of detection was determined with a quantified SARS-CoV-2 virus and has been evaluated at 2.4×10^2 TCID₅₀/mL.

5. Stability Study

5.1 Objective

The stability study was performed to evaluate the shelf life of NANOREPRO SARS-CoV-2 Antigen Rapid Test.

5.2 Parameter

Three kit lots, from three independent productions lots were tested. The finished products of NANOREPRO SARS-CoV-2 Antigen Rapid Test were stored for real time storage stability study in (2-8°C) & (18-30°C), and perform an assessment test each month. For accelerated stability study the finished products were stored at 37°C for 23 days.

Testing procedure and result interpretation were performed according to instruction for use.

5.3 Results

Item Lot	Storage Time	Appearance	Positive Coincidence rate	Negative Coincidence rate
20200509	1 day	PASS	8/8	20/20
	4 day	PASS	8/8	20/20
	7 day	PASS	8/8	20/20
	11 day	PASS	8/8	20/20
	14 day	PASS	8/8	20/20
	17 day	PASS	8/8	20/20
	20 day	PASS	8/8	20/20
20200515	23 day	PASS	8/8	20/20
	1 day	PASS	8/8	20/20
	4 day	PASS	8/8	20/20
	7 day	PASS	8/8	20/20
	11 day	PASS	8/8	20/20
	14 day	PASS	8/8	20/20
	17 day	PASS	8/8	20/20
20200522	20 day	PASS	8/8	20/20
	23 day	PASS	8/8	20/20
	1 day	PASS	8/8	20/20
	4 day	PASS	8/8	20/20
	7 day	PASS	8/8	20/20
	11 day	PASS	8/8	20/20
	14 day	PASS	8/8	20/20
17 day	PASS	8/8	20/20	
20 day	PASS	8/8	20/20	
23 day	PASS	8/8	20/20	

5.4 Conclusion

Through real-time stability and accelerated stability study, it is concluded that the shelf-life of the test is 12 months.

6. Repeatability Study

6.1 Objective

To validate the repeatability of NANOREPRO SARS-CoV-2 Antigen Rapid Test.

6.2 Parameter

Three batches of kits were tested with repetitive reference. The test was repeated 10 times for each batch. The quantity of + shows the degree of color to record the test results.

6.3 Results

LOT Times	20200509		20200515		20200522	
	R1	R2	R1	R2	R1	R2
1	++	+	++	+	++	+
2	++	+	++	+	++	+
3	++	+	++	+	++	+
4	++	+	++	+	++	+
5	++	+	++	+	++	+
6	++	+	++	+	++	+
7	++	+	++	+	++	+
8	++	+	++	+	++	+
9	++	+	++	+	++	+
10	++	+	++	+	++	+

6.4 Conclusion

The test results were positive, and there is no obvious difference in color.

SARS-CoV-2 Antigen Rapid Test

NanoRepro AG

Principle:

SARS-CoV-2 Antigen Detection

Time to Results:

15 Minutes

Specimen:

Nasopharyngeal & Oropharyngeal Swab

User:

Rapid Test for Professional Users

Reference:

Pharmaceutical wholesale, medical specialist trade, laboratory trade, doctors, pharmacies, nursing homes

Product code:

B66000

With the help of the SARS-CoV-2 Antigen Rapid Test an acute infection can be detectable quickly, conveniently and inexpensively. In just 15 minutes is the test result on site. The SARS-CoV-2 viral nucleoprotein antigen is detected through visual interpretation of color in a few simple steps.

PCR Validated Samples

- ✓ **97,0 %** Overall Positive Percent Agreement
- ✓ **98,0 %** Overall Negative Percent Agreement

Quick Results:

- ✓ Interpret results in 15 minutes

Efficient Contact Tracing

- ✓ Immediate results allow for contact tracing & isolation to start sooner

Easy-to-use

1. Collect specimen using included swab

2. Add 23 drops of buffer into the extraction tube
3. Insert swab into extraction tub & mix with solution
4. Add 2 drops of solution into the test's sample well and wait results

Diagnostic Accuracy


SARS-CoV-2 Antigen Rapid Test	Relative Sensitivity	97% (91.5 %-99.4%) *
	Relative Specificity	98% (95.7%-99.3%) *

*95% Konfidenzintervall

Ordering Information


Product Code	Product Name	Contents	Box Size (mm)
B66000	NanoRepro SARS-CoV-2 Antigen Rapid Test	20 Tests/ Kit	225x135x80

1


For each specimen, open the foil pouch just before testing and take the test components, and put it on a clean, level surface. Label the tube with the patient identification. For best results, the assay should be performed within one hour Mix the extraction buffer gently. Add 23 drops to the extraction tube.

2a


Insert the swab into the nostril parallel to the palate. Rotating against the nasal wall (to ensure swab contains cells as well as mucus).

2b

Insert the swab completely from the mouth into the throat, centering on the res part of the throat wall and maxillary tonsils and rub the bilateral throat tonsils and throat wall moderately. Avoid touching the tongue and remove the swab.

3


Insert the swab into the extraction tube. Mix well and squeeze the swab 10-15 times by compressing the walls of the tube against the swab.

4

Roll the swab head against the inner wall of the tubes as you remove it. Try to release as much liquid as possible. Dispose of the used swab in accordance with your biohazard waste disposal protocol.

5


Insert nozzle into sample extraction tube. Invert the tube and add 2 drops of solution into the sample well by the gently squeezing the tube.

6

Read the results at 15-20 minutes.

7


Positive

Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T).

Negative

Only one colored band appears, in the control region (C). No apparent colored band appears in the test region (T).

Invalid

Control band fails to appear. Results from any test which has not produced a control band at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.